


ONO-AE1-329

Chemical Properties

CAS No. : 253350-43-5
 Formula: C₂₂H₃₀O₆S₂
 Molecular Weight: 454.6
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	ONO-AE1-329 is a novel agonist of the prostaglandin PGE2 receptor EP4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1997 mL	10.9987 mL	21.9974 mL
5 mM	0.4399 mL	2.1997 mL	4.3995 mL
10 mM	0.220 mL	1.0999 mL	2.1997 mL
50 mM	0.044 mL	0.220 mL	0.4399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481