


## TIPP

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 146369-65-5   |
| Formula:          | C37H38N4O6  |
| Molecular Weight: | 634.72  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|               |  |
|---------------|--|
| Description   | TIPP is an agent of delta opioid antagonist. |
| Targets(IC50) | Others,Opioid Receptor                       |

## Solubility Information

|            |   |
|------------|---|
| Solubility | DMSO: Soluble,<br>(< 1 mg/ml refers to the product slightly soluble or insoluble) |
|------------|---|

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg      |
|-------|-----------|-----------|-----------|
| 1 mM  | 1.5755 mL | 7.8775 mL | 15.755 mL |
| 5 mM  | 0.3151 mL | 1.5755 mL | 3.151 mL  |
| 10 mM | 0.1575 mL | 0.7877 mL | 1.5755 mL |
| 50 mM | 0.0315 mL | 0.1575 mL | 0.3151 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

### Reference

Lian Q, Jiang W, Cheng Y, Cao H, Liu M, Wang J, Li Y, Song X, Wang F. A novel pentapeptide originated from calf thymus named TIPP shows an inhibitory effect on lung allergic inflammation. *Int Immunopharmacol.* 2015 Feb;24(2):256-266. doi: 10.1016/j.intimp.2014.12.019. Epub 2014 Dec 19. PubMed PMID: 25533504.

Lian Q, Cheng Y, Zhong C, Wang F. Inhibition of the IgE-mediated activation of RBL-2H3 cells by TIPP, a novel thymic immunosuppressive pentapeptide. *Int J Mol Sci.* 2015 Jan 20;16(1):2252-68. doi: 10.3390/ijms16012252. PubMed PMID: 25608657; PubMed Central PMCID: PMC4307361.

Koning GG, Adang EM, Stalmeier PF, Keus F, Vriens PW, van Laarhoven CJ. TIPP and Lichtenstein modalities for inguinal hernia repair: a cost minimisation analysis alongside a randomised trial. *Eur J Health Econ.* 2013 Dec;14(6):1027-34. doi: 10.1007/s10198-012-0453-0. Epub 2012 Dec 28. PubMed PMID: 23271350.

Nguyen NP, Mirarab S, Liu B, Pop M, Warnow T. TIPP: taxonomic identification and phylogenetic profiling. *Bioinformatics.* 2014 Dec 15;30(24):3548-55. doi: 10.1093/bioinformatics/btu721. Epub 2014 Oct 29. PubMed PMID: 25359891; PubMed Central PMCID: PMC4253836.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481