


## Erysoitrine

## Chemical Properties

CAS No. :	27740-43-8
Formula:	C <sub>19</sub> H <sub>23</sub> NO <sub>3</sub>
Molecular Weight:	313.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Erysoitrine, shows antibacterial activities, isolated from seed pods of Erythrina latissima.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1909 mL	15.9546 mL	31.9091 mL
5 mM	0.6382 mL	3.1909 mL	6.3818 mL
10 mM	0.3191 mL	1.5955 mL	3.1909 mL
50 mM	0.0638 mL	0.3191 mL	0.6382 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

de Ávila JM, et al. Antimicrobial Evaluation of Erythrinan Alkaloids from Erythrina cristagalli L. Med Chem. 2018;14 (8):784-790.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481