


5-Methoxy-2-tetralone

Chemical Properties

CAS No. : 32940-15-1
 Formula: C₁₁H₁₂O₂
 Molecular Weight: 176.215
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Methoxy-2-tetralone, with CAS No. 32940-15-1, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 5-Methoxy-2-tetralone provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6747 mL	28.3736 mL	56.7472 mL
5 mM	1.1349 mL	5.6747 mL	11.3494 mL
10 mM	0.5675 mL	2.8374 mL	5.6747 mL
50 mM	0.1135 mL	0.5675 mL	1.1349 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481