


## MRS 2957 triethylammonium salt

## Chemical Properties

CAS No. :	1228271-30-4
Formula:	C <sub>19</sub> H <sub>28</sub> N <sub>5</sub> O <sub>20</sub> P <sub>3</sub>
Molecular Weight:	739.37
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	MRS 2957 triethylammonium salt is a P2Y6 agonist.
Targets(IC50)	Others,P2Y Receptor

## Solubility Information

Solubility	H <sub>2</sub> O: <10 mM,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3525 mL	6.7625 mL	13.525 mL
5 mM	0.2705 mL	1.3525 mL	2.705 mL
10 mM	0.1353 mL	0.6763 mL	1.3525 mL
50 mM	0.0271 mL	0.1353 mL	0.2705 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481