

Boc-NH-PEG1-CH₂COOH

Chemical Properties

CAS No. : 142929-49-5


Formula: C₉H₁₇NO₅

Molecular Weight: 219.24

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boc-NH-PEG1-CH ₂ COOH is a PEG-based linker for PROTACs that connects two essential ligands, facilitating selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binding to an E3 ubiquitin ligase and the other to the target protein. By leveraging the intracellular ubiquitin-proteasome system, PROTACs enable the selective degradation of target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5612 mL	22.8061 mL	45.6121 mL
5 mM	0.9122 mL	4.5612 mL	9.1224 mL
10 mM	0.4561 mL	2.2806 mL	4.5612 mL
50 mM	0.0912 mL	0.4561 mL	0.9122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481