


L-Lysyl-L-lysineamide

Chemical Properties

CAS No. :	134276-59-8
Formula:	C ₁₂ H ₂₇ N ₅ O ₂
Molecular Weight:	273.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-Lysyl-L-lysineamide is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6579 mL	18.2896 mL	36.5791 mL
5 mM	0.7316 mL	3.6579 mL	7.3158 mL
10 mM	0.3658 mL	1.829 mL	3.6579 mL
50 mM	0.0732 mL	0.3658 mL	0.7316 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Drugs and Lactation Database (LactMed) [Internet]. Bethesda (MD): National Library of Medicine (US); 2006-. Available from <http://www.ncbi.nlm.nih.gov/books/NBK500975/> PubMed PMID: 30000034.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481