


Azido-PEG4-Boc

Chemical Properties

CAS No. :	581066-04-8
Formula:	C15H29N3O6
Molecular Weight:	347.41
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Azido-PEG4-Boc is a PEG-based linker for PROTACs that joins two essential ligands, pivotal for forming PROTAC molecules, enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binds to an E3 ubiquitin ligase, and the other targets a specific protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8784 mL	14.3922 mL	28.7844 mL
5 mM	0.5757 mL	2.8784 mL	5.7569 mL
10 mM	0.2878 mL	1.4392 mL	2.8784 mL
50 mM	0.0576 mL	0.2878 mL	0.5757 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481