

Compound 0080-0014

## Chemical Properties


CAS No. :

Formula: C<sub>28</sub>H<sub>32</sub>O<sub>14</sub>

Molecular Weight: 592.55

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0080-0014 is a useful organic compound for research related to life sciences and the catalog number is T131705.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6876 mL	8.4381 mL	16.8762 mL
5 mM	0.3375 mL	1.6876 mL	3.3752 mL
10 mM	0.1688 mL	0.8438 mL	1.6876 mL
50 mM	0.0338 mL	0.1688 mL	0.3375 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481