


Amipizone

Chemical Properties

CAS No. : 69635-63-8
 Formula: C₁₄H₁₆ClN₃O₂
 Molecular Weight: 293.75
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Amipizone is a phosphodiesterase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4043 mL	17.0213 mL	34.0426 mL
5 mM	0.6809 mL	3.4043 mL	6.8085 mL
10 mM	0.3404 mL	1.7021 mL	3.4043 mL
50 mM	0.0681 mL	0.3404 mL	0.6809 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481