

Antidepressant agent 2

Chemical Properties

CAS No. :

Formula: C₂₁H₂₂ClFN₂O₃S

Molecular Weight: 436.93

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Antidepressant agent 2 showed significant antidepressant effects with a MED value of 0.1 mg/kg.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2887 mL	11.4435 mL	22.887 mL
5 mM	0.4577 mL	2.2887 mL	4.5774 mL
10 mM	0.2289 mL	1.1443 mL	2.2887 mL
50 mM	0.0458 mL	0.2289 mL	0.4577 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481