

## 2-(2-hydroxy-3-methylbutanamido)-3-methylbutanoic acid

## Chemical Properties

CAS No. :

Formula: C<sub>10</sub>H<sub>19</sub>NO<sub>4</sub>

Molecular Weight: 217.27

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	2-(2-hydroxy-3-methylbutanamido)-3-methylbutanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6026 mL	23.0128 mL	46.0257 mL
5 mM	0.9205 mL	4.6026 mL	9.2051 mL
10 mM	0.4603 mL	2.3013 mL	4.6026 mL
50 mM	0.0921 mL	0.4603 mL	0.9205 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481