

## Ginsenoyne K

## Chemical Properties

CAS No. :	141947-42-4
Formula:	C17H24O3
Molecular Weight:	276.37
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Ginsenoyne K is a useful organic compound for research related to life sciences and the catalog number is T125902.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6183 mL	18.0917 mL	36.1834 mL
5 mM	0.7237 mL	3.6183 mL	7.2367 mL
10 mM	0.3618 mL	1.8092 mL	3.6183 mL
50 mM	0.0724 mL	0.3618 mL	0.7237 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481