

AX15892

Chemical Properties


CAS No. :

Formula: C32H40N8O3

Molecular Weight: 584.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	AX15892 is a potent and selective inhibitor of ERK5.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7102 mL	8.551 mL	17.1019 mL
5 mM	0.342 mL	1.7102 mL	3.4204 mL
10 mM	0.171 mL	0.8551 mL	1.7102 mL
50 mM	0.0342 mL	0.171 mL	0.342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481