


## APTO-253 isomer

## Chemical Properties

CAS No. :	1422826-80-9
Formula:	C <sub>22</sub> H <sub>14</sub> FN <sub>5</sub>
Molecular Weight:	367.38
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	APTO-253 isomer is an isomer of APTO-253 with the Fluoro-atom in 6-position. APTO-253 also known as LOR-253, LT-253, is a small molecule inhibitor of human metal-regulatory transcription factor 1 (MTF-1) with potential antitumor activity. MTF-1 inhibitor LOR-253 inhibits MTF-1 activity and thereby induces the expression of MTF-1 dependent tumor suppressor factor Kruppel like factor 4 (KLF4). This subsequently leads to the downregulation of cyclin D1, blocking cell cycle progression and proliferation. This agent also causes decreased expression of genes involved in tumor hypoxia and angiogenesis.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.722 mL	13.6099 mL	27.2198 mL
5 mM	0.5444 mL	2.722 mL	5.444 mL
10 mM	0.2722 mL	1.361 mL	2.722 mL
50 mM	0.0544 mL	0.2722 mL	0.5444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481