

## (8E,12Z)-10-hydroxyoctadeca-8,12-dienoic acid

## Chemical Properties

CAS No. :

Formula: C<sub>18</sub>H<sub>32</sub>O<sub>3</sub>

Molecular Weight: 296.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(8E,12Z)-10-hydroxyoctadeca-8,12-dienoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3733 mL	16.8663 mL	33.7325 mL
5 mM	0.6747 mL	3.3733 mL	6.7465 mL
10 mM	0.3373 mL	1.6866 mL	3.3733 mL
50 mM	0.0675 mL	0.3373 mL	0.6747 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481