

3-(2,4-dihydroxypentyl)-6,8-dihydroxy-1H-isochromen-1-one

Chemical Properties


CAS No. :

Formula: C₁₄H₁₆O₆

Molecular Weight: 280.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

3-(2,4-dihydroxypentyl)-6,8-dihydroxy-1H-isochromen-1-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5679 mL	17.8393 mL	35.6786 mL
5 mM	0.7136 mL	3.5679 mL	7.1357 mL
10 mM	0.3568 mL	1.7839 mL	3.5679 mL
50 mM	0.0714 mL	0.3568 mL	0.7136 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481