


cis-Hexahydrophthalimide

Chemical Properties

CAS No. : 7506-66-3
 Formula: C₈H₁₁NO₂
 Molecular Weight: 147.13
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	cis-Hexahydrophthalimide, with CAS No. 7506-66-3, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. cis-Hexahydrophthalimide provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7967 mL	33.9836 mL	67.9671 mL
5 mM	1.3593 mL	6.7967 mL	13.5934 mL
10 mM	0.6797 mL	3.3984 mL	6.7967 mL
50 mM	0.1359 mL	0.6797 mL	1.3593 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481