


2MeSAMP

Chemical Properties

CAS No. :	81921-45-1
Formula:	C ₁₁ H ₁₆ N ₅ NaO ₇ PS
Molecular Weight:	416.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2MeSAMP is a P2Y ₁₂ antagonist that works by inhibiting platelet activation through a P2Y ₁₂ /G(i)-dependent mechanism.
Targets(IC50)	Others,P2Y Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4021 mL	12.0106 mL	24.0211 mL
5 mM	0.4804 mL	2.4021 mL	4.8042 mL
10 mM	0.2402 mL	1.2011 mL	2.4021 mL
50 mM	0.048 mL	0.2402 mL	0.4804 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Tourdot BE, Conaway S, Niisuke K, Edelstein LC, Bray PF, Holinstat M. Mechanism of race-dependent platelet activation through the protease-activated receptor-4 and Gq signaling axis. *Arterioscler Thromb Vasc Biol.* 2014 Dec;34(12):2644-50. doi: 10.1161/ATVBAHA.114.304249. Epub 2014 Oct 2. PubMed PMID: 25278289; PubMed Central PMCID: PMC4239175.

Li R, Diamond SL. Detection of platelet sensitivity to inhibitors of COX-1, P2Y₁, and P2Y₁₂ using a whole blood microfluidic flow assay. *Thromb Res.* 2014 Feb;133(2):203-10. doi: 10.1016/j.thromres.2013.10.043. Epub 2013 Nov 6. PubMed PMID: 24365044; PubMed Central PMCID: PMC3932538.

Liñán-Rico A, Wunderlich JE, Grants IS, Frankel WL, Xue J, Williams KC, Harzman AE, Enneking JT, Cooke HJ, Christofi FL. Purinergic autocrine regulation of mechanosensitivity and serotonin release in a human EC model: ATP-gated P2X3 channels in EC are downregulated in ulcerative colitis. *Inflamm Bowel Dis.* 2013 Oct;19(11):2366-79. doi: 10.1097/MIB.0b013e31829ecf4d. PubMed PMID: 23917247; PubMed Central PMCID: PMC4037929.

Herfindal L, Nygaard G, Kopperud R, Krakstad C, Døskeland SO, Selheim F. Off-target effect of the Epac agonist 8-pCPT-2'-O-Me-cAMP on P2Y₁₂ receptors in blood platelets. *Biochem Biophys Res Commun.* 2013 Aug 9;437(4):603-8. doi: 10.1016/j.bbrc.2013.07.007. Epub 2013 Jul 12. PubMed PMID: 23850619.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481