


2'-O-Methylisoliquiritigenin

Chemical Properties

CAS No. :	51828-10-5
Formula:	C ₁₆ H ₁₄ O ₄
Molecular Weight:	270.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-O-Methylisoliquiritigenin, a compound synthesized by enzymes specifically induced in NR.2'-O-Methylisoliquiritigenin, isolated from the Arachis species, up-regulates 5-HT, NE, DA and GABA pathways, but does not put a very significant effect on ne NE pathway.
Targets(IC50)	5-HT Receptor,GABA Receptor,Dopamine Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6999 mL	18.4993 mL	36.9987 mL
5 mM	0.740 mL	3.6999 mL	7.3997 mL
10 mM	0.370 mL	1.8499 mL	3.6999 mL
50 mM	0.074 mL	0.370 mL	0.740 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lei Deng, et al. Sedative-hypnotic and Anxiolytic Effects and the Mechanism of Action of Aqueous Extracts of Peanut Stems and Leaves in Mice. J Sci Food Agric. 2018 Oct;98(13):4885-4894.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481