

4-{2,3-dihydroxy-5-methyl-7-oxabicyclo[2.2.1]hept-5-en-2-yl}-3-methylbutanoic acid

Chemical Properties


CAS No. :

Formula: C₁₂H₁₈O₅

Molecular Weight: 242.27

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-{2,3-dihydroxy-5-methyl-7-oxabicyclo[2.2.1]hept-5-en-2-yl}-3-methylbutanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1276 mL	20.6381 mL	41.2763 mL
5 mM	0.8255 mL	4.1276 mL	8.2553 mL
10 mM	0.4128 mL	2.0638 mL	4.1276 mL
50 mM	0.0826 mL	0.4128 mL	0.8255 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481