

DBCO-PEG2-NH-Boc

Chemical Properties

CAS No. : 2377004-08-3


Formula: C30H37N3O6

Molecular Weight: 535.63

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	DBCO-PEG2-NH-Boc is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs comprise two distinct ligands joined by a linker: one ligand binds to an E3 ubiquitin ligase, and the other targets a specific protein. These compounds utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.867 mL	9.3348 mL	18.6696 mL
5 mM	0.3734 mL	1.867 mL	3.7339 mL
10 mM	0.1867 mL	0.9335 mL	1.867 mL
50 mM	0.0373 mL	0.1867 mL	0.3734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481