


Bradykinin potentiator-5

Chemical Properties

CAS No. :	30505-63-6
Formula:	C30H41N7O7
Molecular Weight:	611.69
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bradykinin potentiator-5 is a peptide inflammatory mediator, causes blood vessels to dilate (enlarge), and therefore causes blood pressure to fall.
Targets(IC50)	RAAS,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6348 mL	8.1741 mL	16.3482 mL
5 mM	0.327 mL	1.6348 mL	3.2696 mL
10 mM	0.1635 mL	0.8174 mL	1.6348 mL
50 mM	0.0327 mL	0.1635 mL	0.327 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Balaguer JM, Yu C, Byrne JG, Ball SK, Petracek MR, Brown NJ, Pretorius M. Contribution of endogenous bradykinin to fibrinolysis, inflammation, and blood product transfusion following cardiac surgery: a randomized clinical trial. *Clin Pharmacol Ther.* 2013 Apr;93(4):326-34. doi: 10.1038/clpt.2012.249. Epub 2012 Dec 24. PubMed PMID: 23361105; PubMed Central PMCID: PMC4031681.

Jiang Y, Xi X, Ge L, Yang N, Hou X, Ma J, Ma C, Wu Y, Guo X, Li R, Zhou M, Wang L, Chen T, Shaw C. Bradykinin-related peptides (BRPs) from skin secretions of three genera of phyllomedusine leaf frogs and their comparative pharmacological effects on mammalian smooth muscles. *Peptides.* 2014 Feb;52:122-33. doi: 10.1016/j.peptides.2013.12.013. Epub 2014 Jan 3. PubMed PMID: 24394432.

Shi D, Luo Y, Du Q, Wang L, Zhou M, Ma J, Li R, Chen T, Shaw C. A novel bradykinin-related dodecapeptide (RVALPPGFPLR) from the skin secretion of the Fujian large-headed frog (*Limnonectes fujianensis*) exhibiting unusual structural and functional features. *Toxins (Basel).* 2014 Sep 29;6(10):2886-98. doi: 10.3390/toxins6102886. PubMed PMID: 25268979; PubMed Central PMCID: PMC4210874.

Uzawa A, Mori M, Taniguchi J, Kuwabara S. Modulation of the kallikrein/kinin system by the angiotensin-converting enzyme inhibitor alleviates experimental autoimmune encephalomyelitis. *Clin Exp Immunol.* 2014 Nov;178(2):245-52. doi: 10.1111/cei.12413. PubMed PMID: 24996009; PubMed Central PMCID: PMC4233374.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481