


## Tenidap Sodium

## Chemical Properties

CAS No. :	119784-94-0
Formula:	C <sub>14</sub> H <sub>8</sub> ClN <sub>2</sub> NaO <sub>3</sub> S
Molecular Weight:	342.72
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tenidap Sodium is a COX/5-LOX inhibitor and cytokine-modulating anti-inflammatory agent.
Targets(IC50)	Others,COX

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9178 mL	14.5892 mL	29.1783 mL
5 mM	0.5836 mL	2.9178 mL	5.8357 mL
10 mM	0.2918 mL	1.4589 mL	2.9178 mL
50 mM	0.0584 mL	0.2918 mL	0.5836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481