

Dapiglutide

Chemical Properties

CAS No. : 2296814-85-0

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Dapiglutide

Biological Description

Description	Dapiglutide (ZP7570) is a long-acting dual agonist of both glucagon-like peptide-1 receptor (GLP-1R) and glucagon-like peptide-2 receptor (GLP-2R), with potential applications in short bowel syndrome (SBS) research.
Targets(IC50)	Glucagon Receptor

Reference

Zealand Pharma doses first subject in Phase 1 clinical trial with potential next generation dual-acting peptide therapeutic for short bowel syndrome. 26 June 2019.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481