


Benzyl p-coumarate

Chemical Properties

CAS No. :	61844-62-0
Formula:	C16H14O3
Molecular Weight:	254.285
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzyl p-coumarate is a natural product from poplar resins.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9325 mL	19.6626 mL	39.3252 mL
5 mM	0.7865 mL	3.9325 mL	7.865 mL
10 mM	0.3933 mL	1.9663 mL	3.9325 mL
50 mM	0.0787 mL	0.3933 mL	0.7865 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481