


Eglandulosine

Chemical Properties

CAS No. :	53508-33-1
Formula:	C ₂₁ H ₂₄ N ₂ O ₃
Molecular Weight:	352.434
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Eglandulosine is a useful organic compound for research related to life sciences. The catalog number is T124699 and the CAS number is 53508-33-1.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8374 mL	14.1872 mL	28.3744 mL
5 mM	0.5675 mL	2.8374 mL	5.6749 mL
10 mM	0.2837 mL	1.4187 mL	2.8374 mL
50 mM	0.0567 mL	0.2837 mL	0.5675 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481