


## FGF/PDGF/VEGF RTK Inhibitor

## Chemical Properties

CAS No. :	144335-37-5
Formula:	C <sub>16</sub> H <sub>11</sub> N <sub>3</sub> O <sub>2</sub>
Molecular Weight:	277.28
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	FGF/PDGF/VEGF RTK Inhibitor is a potent, reversible, atp-competitive inhibitor against pdgfrβ, fgfr-1, and vegfr-2 (ic50 = 20, 90, and 240 nm, respectively) and effectively suppresses vegf-stimulated proliferation of hmvecs (ec50 = 420 nm).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6065 mL	18.0323 mL	36.0646 mL
5 mM	0.7213 mL	3.6065 mL	7.2129 mL
10 mM	0.3606 mL	1.8032 mL	3.6065 mL
50 mM	0.0721 mL	0.3606 mL	0.7213 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481