


5-HT3 antagonist 1

Chemical Properties

CAS No. :	129294-09-3
Formula:	C ₂₂ H ₂₇ N ₅ O
Molecular Weight:	377.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-HT3 antagonist 1 is a potent and selective inhibitor of the 5-HT3 receptor.
Targets(IC50)	Others,5-HT Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6491 mL	13.2457 mL	26.4915 mL
5 mM	0.5298 mL	2.6491 mL	5.2983 mL
10 mM	0.2649 mL	1.3246 mL	2.6491 mL
50 mM	0.053 mL	0.2649 mL	0.5298 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

INDAZOLE-3-CARBOXYLIC ACID DERIVATIVES. US5017573

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481