


## 3,4-Dihydroxy-2-O-methylanigorufone

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 1392307-42-4  |
| Formula:          | C <sub>20</sub> H <sub>14</sub> O <sub>4</sub>  |
| Molecular Weight: | 318.32  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|               |  |
|---------------|--|
| Description   | 3,4-Dihydroxy-2-O-methylanigorufone is a natural product for research related to life sciences. The catalog number is TN2877 and the CAS number is 1392307-42-4. |
| Targets(IC50) | Others   |

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 3.1415 mL | 15.7075 mL | 31.4149 mL |
| 5 mM  | 0.6283 mL | 3.1415 mL  | 6.283 mL   |
| 10 mM | 0.3141 mL | 1.5707 mL  | 3.1415 mL  |
| 50 mM | 0.0628 mL | 0.3141 mL  | 0.6283 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481