


## (S)-GNA-C(Bz)-phosphoramidite

## Chemical Properties

CAS No. :	885322-06-5
Formula:	C44H50N5O7P
Molecular Weight:	791.87
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-GNA-C(Bz)-phosphoramidite is a Nucleoside Derivative - Acyclic nucleoside; Nucleoside Phosphoramidite.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2628 mL	6.3142 mL	12.6283 mL
5 mM	0.2526 mL	1.2628 mL	2.5257 mL
10 mM	0.1263 mL	0.6314 mL	1.2628 mL
50 mM	0.0253 mL	0.1263 mL	0.2526 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481