


Moexipril ethyl ester

Chemical Properties

CAS No. :	103733-37-5
Formula:	C ₂₉ H ₃₈ N ₂ O ₇
Molecular Weight:	526.62
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Moexipril ethyl ester is a long-acting, non-sulphydryl angiotensine-converting enzyme inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8989 mL	9.4945 mL	18.989 mL
5 mM	0.3798 mL	1.8989 mL	3.7978 mL
10 mM	0.1899 mL	0.9495 mL	1.8989 mL
50 mM	0.038 mL	0.1899 mL	0.3798 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481