


## 2-(1,2,3,4-Tetrahydroisoquinolin-1-yl)acetic acid hydrochloride

## Chemical Properties

CAS No. :	53014-68-9
Formula:	C <sub>11</sub> H <sub>14</sub> ClNO <sub>2</sub>
Molecular Weight:	227.687
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-(1,2,3,4-Tetrahydroisoquinolin-1-yl)acetic acid hydrochloride, with CAS No. 53014-68-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-(1,2,3,4-Tetrahydroisoquinolin-1-yl)acetic acid hydrochloride provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3919 mL	21.9597 mL	43.9194 mL
5 mM	0.8784 mL	4.3919 mL	8.7839 mL
10 mM	0.4392 mL	2.196 mL	4.3919 mL
50 mM	0.0878 mL	0.4392 mL	0.8784 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481