


## Arfendazam

## Chemical Properties

CAS No. :	37669-57-1
Formula:	C <sub>18</sub> H <sub>17</sub> ClN <sub>2</sub> O <sub>3</sub>
Molecular Weight:	344.79
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Arfendazam is a Tranquilizer.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9003 mL	14.5016 mL	29.0032 mL
5 mM	0.5801 mL	2.9003 mL	5.8006 mL
10 mM	0.290 mL	1.4502 mL	2.9003 mL
50 mM	0.058 mL	0.290 mL	0.5801 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481