


5-[3-[(Trifluoroacetyl)amino]propyl]uridine

Chemical Properties

CAS No. :	327623-80-3
Formula:	C ₁₄ H ₁₈ F ₃ N ₃ O ₇
Molecular Weight:	397.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-[3-[(Trifluoroacetyl)amino]propyl]uridine is a nucleoside derivative, specifically a 5-modified pyrimidine nucleoside, and is categorized as an amino nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.517 mL	12.5849 mL	25.1699 mL
5 mM	0.5034 mL	2.517 mL	5.034 mL
10 mM	0.2517 mL	1.2585 mL	2.517 mL
50 mM	0.0503 mL	0.2517 mL	0.5034 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481