

## Poly(4-vinylphenol)

### Chemical Properties

CAS No. : 24979-70-2

Formula: (C<sub>8</sub>H<sub>8</sub>O)<sub>x</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Poly(4-vinylphenol) is an endogenous metabolite.
Targets(IC50)	Others,Endogenous Metabolite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481