


Effusanin E

Chemical Properties

CAS No. :	76470-15-0
Formula:	C ₂₀ H ₂₈ O ₆
Molecular Weight:	364.438
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Effusanin E shows antibacterial and anti-cancer activities, it significantly inhibits cell proliferation and induces apoptosis in NPC cells by suppressing p50/p65 proteins to down-regulate COX-2 expression .
Targets(IC50)	NF-κB,Caspase,Antifection,COX,PARP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7439 mL	13.7197 mL	27.4394 mL
5 mM	0.5488 mL	2.7439 mL	5.4879 mL
10 mM	0.2744 mL	1.372 mL	2.7439 mL
50 mM	0.0549 mL	0.2744 mL	0.5488 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481