


Boc-D-Phe(3,4-DiF)-OH

Chemical Properties

CAS No. :	205445-51-8
Formula:	C ₁₄ H ₁₇ F ₂ N ₁ O ₄
Molecular Weight:	301.29
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Boc-D-Phe(3,4-DiF)-OH, an amino acid derivative, has extensive applications in life science-related research.
Targets(IC50)	Amino Acids and Derivatives

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3191 mL	16.5953 mL	33.1906 mL
5 mM	0.6638 mL	3.3191 mL	6.6381 mL
10 mM	0.3319 mL	1.6595 mL	3.3191 mL
50 mM	0.0664 mL	0.3319 mL	0.6638 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481