

## Compound T129174

## Chemical Properties


CAS No. :

Formula: C42H64O15

Molecular Weight: 808.96

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3S,4S,5R,6R)-6-(((3S,6aR,6bS,8aR,9R,14bS)-14b-formyl-9-hydroxy-4,4,6a,6b,8a,11,11-heptamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,14,14a,14b-icosahydricen-3-yl]oxy)-5-(((2R,3R,4S,5S,6S)-6-carboxy-3,4,5-trihydroxyoxan-2-yl]oxy)-3,4-dihydroxyoxane-2-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2362 mL	6.1808 mL	12.3616 mL
5 mM	0.2472 mL	1.2362 mL	2.4723 mL
10 mM	0.1236 mL	0.6181 mL	1.2362 mL
50 mM	0.0247 mL	0.1236 mL	0.2472 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481