

## Azido-PEG3-methyl ester

## Chemical Properties

CAS No. : 2086689-07-6


Formula: C10H19N3O5

Molecular Weight: 261.28

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Azido-PEG3-methyl ester is a polyethylene glycol (PEG)-derived linker utilized in the synthesis of proteolysis-targeting chimeras (PROTACs)[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands connected by a linker: one binds to an E3 ubiquitin ligase and the other to the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8273 mL	19.1366 mL	38.2731 mL
5 mM	0.7655 mL	3.8273 mL	7.6546 mL
10 mM	0.3827 mL	1.9137 mL	3.8273 mL
50 mM	0.0765 mL	0.3827 mL	0.7655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481