

Compound 0449-0140

## Chemical Properties


CAS No. :

Formula: C<sub>35</sub>H<sub>54</sub>O<sub>13</sub>

Molecular Weight: 682.804

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound 0449-0140 is a useful organic compound for research related to life sciences and the catalog number is T131645.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4646 mL	7.3228 mL	14.6456 mL
5 mM	0.2929 mL	1.4646 mL	2.9291 mL
10 mM	0.1465 mL	0.7323 mL	1.4646 mL
50 mM	0.0293 mL	0.1465 mL	0.2929 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481