


L 652117

Chemical Properties

CAS No. : 127020-81-9
 Formula: C₁₄H₁₆N₂O₇S
 Molecular Weight: 356.35
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L 652117 is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8062 mL	14.0311 mL	28.0623 mL
5 mM	0.5612 mL	2.8062 mL	5.6125 mL
10 mM	0.2806 mL	1.4031 mL	2.8062 mL
50 mM	0.0561 mL	0.2806 mL	0.5612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Wilmouth RC, Westwood NJ, Anderson K, Brownlee W, Claridge TD, Clifton IJ, Pritchard GJ, Aplin RT, Schofield CJ. Inhibition of elastase by N-sulfonylaryl beta-lactams: anatomy of a stable acyl-enzyme complex. *Biochemistry*. 1998 Dec 15;37(50):17506-13. PubMed PMID: 9860865.
- Knight WB, Swiderek KM, Sakuma T, Calaycay J, Shively JE, Lee TD, Covey TR, Shushan B, Green BG, Chabin R, et al. Electrospray ionization mass spectrometry as a mechanistic tool: mass of human leucocyte elastase and a beta-lactam-derived E-I complex. *Biochemistry*. 1993 Mar 2;32(8):2031-5. PubMed PMID: 8095408.
- Knight WB, Chabin R, Green B. Inhibition of human serine proteases by substituted 2-azetidinones. *Arch Biochem Biophys*. 1992 Aug 1;296(2):704-8. PubMed PMID: 1632655.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481