


Deoxymikanolide

Chemical Properties

CAS No. : 23753-57-3
 Formula: C₁₅H₁₆O₅
 Molecular Weight: 276.288
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Deoxymikanolide is a natural product for research related to life sciences. The catalog number is TN5393 and the CAS number is 23753-57-3.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6194 mL	18.0969 mL	36.1939 mL
5 mM	0.7239 mL	3.6194 mL	7.2388 mL
10 mM	0.3619 mL	1.8097 mL	3.6194 mL
50 mM	0.0724 mL	0.3619 mL	0.7239 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481