

PEG6-(CH₂CO₂H)₂

Chemical Properties

CAS No. : 77855-76-6


Formula: C₁₄H₂₆O₁₀

Molecular Weight: 354.35

Keep away from direct sunlight

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PEG6-(CH ₂ CO ₂ H) ₂ is a symmetric PEG PROTAC linker.
Targets(IC ₅₀)	E3 Ligase Ligand-Linker Conjugates, PROTAC Linker

Solubility Information

Solubility	DMSO: 10 mM, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8221 mL	14.1103 mL	28.2207 mL
5 mM	0.5644 mL	2.8221 mL	5.6441 mL
10 mM	0.2822 mL	1.411 mL	2.8221 mL
50 mM	0.0564 mL	0.2822 mL	0.5644 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Maniaci C, et al. Homo-PROTACs: bivalent small-molecule dimerizers of the VHL E3 ubiquitin ligase to induce self-degradation. Nat Commun. 2017 Oct 10;8(1):830.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481