

Siwanine A

Chemical Properties


CAS No. :

Formula: C₂₇H₃₉N₀O₈

Molecular Weight: 505.608

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Siwanine A is a useful organic compound for research related to life sciences and the catalog number is T125700.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9778 mL	9.889 mL	19.7781 mL
5 mM	0.3956 mL	1.9778 mL	3.9556 mL
10 mM	0.1978 mL	0.9889 mL	1.9778 mL
50 mM	0.0396 mL	0.1978 mL	0.3956 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481