


## Xylostinine

## Chemical Properties

CAS No. :	50474-67-4
Formula:	C17H34N4O10
Molecular Weight:	454.47
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Xylostinine is a biochemical.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2004 mL	11.0018 mL	22.0037 mL
5 mM	0.4401 mL	2.2004 mL	4.4007 mL
10 mM	0.220 mL	1.1002 mL	2.2004 mL
50 mM	0.044 mL	0.220 mL	0.4401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481