


TrxR1-IN-B19

## Chemical Properties

CAS No. : 170950-29-5  
 Formula: C<sub>21</sub>H<sub>22</sub>O<sub>5</sub>  
 Molecular Weight: 354.4  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	TrxR1-IN-B19 (TrxR1 IN B19), a TrxR1 inhibitor, selectively kills gastric cancer cells by a small molecule via targeting TrxR1 and ROS-mediated ER stress activation.
Targets(IC50)	Apoptosis,Reductase

## Solubility Information

Solubility	DMSO: 65 mg/mL (183.41 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8217 mL	14.1084 mL	28.2167 mL
5 mM	0.5643 mL	2.8217 mL	5.6433 mL
10 mM	0.2822 mL	1.4108 mL	2.8217 mL
50 mM	0.0564 mL	0.2822 mL	0.5643 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481