

Ac-ATS010-KE

## Chemical Properties


CAS No. : 2379342-67-1

Formula: C43H41F5N6O12S

Molecular Weight: 960.88

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ac-ATS010-KE is a novel Selective, irreversible, and Rapid Cell-Permeable Inhibitor of Human Caspase-3.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0407 mL	5.2036 mL	10.4071 mL
5 mM	0.2081 mL	1.0407 mL	2.0814 mL
10 mM	0.1041 mL	0.5204 mL	1.0407 mL
50 mM	0.0208 mL	0.1041 mL	0.2081 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481