


8,8''-Biskoenigine

Chemical Properties

CAS No. :	477890-82-7
Formula:	C38H36N2O6
Molecular Weight:	616.7
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	8,8''-Biskoenigine is a natural product for research related to life sciences. The catalog number is TN6197 and the CAS number is 477890-82-7.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6215 mL	8.1077 mL	16.2153 mL
5 mM	0.3243 mL	1.6215 mL	3.2431 mL
10 mM	0.1622 mL	0.8108 mL	1.6215 mL
50 mM	0.0324 mL	0.1622 mL	0.3243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481