


AX-7396

## Chemical Properties

CAS No. : 824431-12-1  
 Formula: C<sub>15</sub>H<sub>11</sub>Cl<sub>2</sub>N<sub>3</sub>  
 Molecular Weight: 304.17  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	AX-7396 is a protein kinase inhibitor with anti-cytomegaloviral properties.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2876 mL	16.4382 mL	32.8764 mL
5 mM	0.6575 mL	3.2876 mL	6.5753 mL
10 mM	0.3288 mL	1.6438 mL	3.2876 mL
50 mM	0.0658 mL	0.3288 mL	0.6575 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481